




# Gynemed - Microtech Micropipettes










**Safe**.....your success

**Complete**.....your demands

**Easy**.....your administration

<b>Holding pipettes</b>				
<b>Identity colour: Green</b>	<b>Product information</b>	<b>Ordering information</b>	<b>Article number</b>	<b>Box</b>
	<b>Holding micropipettes</b> <ul style="list-style-type: none"> <li>- prepared from borosilicate glass tubing (O.D. 1mm, I.D. 0.75mm)</li> <li>- total length of the pipettes 5.5cm</li> <li>- bending angle 20°-45°</li> <li>- length of tip 1mm,</li> <li>- opening 15-25µm I.D., polished.</li> </ul>	O.D. 80µm, bend 35 degrees, polished opening 15µm	001-80-35	20
		O.D. 100µm, bend 35 degrees, polished opening 15µm	001-100-35	20
		O.D. 120µm, bend 35 degrees, polished opening 15µm	001-120-35	20
		O.D. 80µm, bend 30 degrees, polished opening 15µm	001-80-30	20
		O.D. 100µm, bend 30 degrees, polished opening 15µm	001-100-30	20
		O.D. 120µm, bend 30 degrees, polished opening 15µm	001-120-30	20
		<b>Injection pipettes</b>		
<b>Identity colour: Blue</b>	<b>Product information</b>	<b>Ordering information</b>	<b>Article number</b>	<b>Box</b>
	<b>Injection micropipettes</b> <ul style="list-style-type: none"> <li>- prepared from borosilicate glass tubing (O.D. 1mm, I.D. 0.78mm),</li> <li>- total length of the pipettes 5.5cm,</li> <li>- bending angle 20°-45°</li> <li>- length of tip 0.5mm</li> <li>- bevelled 35°</li> <li>- with or without spike</li> <li>- with I.D. of the tip 4.5-5µm.</li> </ul>	I.D. 4.5-5µm, O.D. 6-7µm, bend 35 degrees, with spike (s)	002-5-35-s	20
		I.D. 4.5-5µm, O.D. 6-7µm, bend 35 degrees, without spike (s)	002-5-35	20
		I.D. 4.5-5µm, O.D. 6-7µm, bend 30 degrees, with spike (s)	002-5-30-s	20
		I.D. 4.5-5µm, O.D. 6-7µm, bend 30 degrees, without spike (s)	002-5-30	20
		I.D. 4.5-5µm, O.D. 6-7µm, bend 20 degrees, with spike (s)	002-5-20-s	20
		I.D. 4.5-5µm, O.D. 6-7µm, bend 20 degrees, without spike (s)	002-5-20	20
		I.D. 4.5-5µm, O.D. 6-7µm, with spike (s)	002-5-s	20
		I.D. 4.5µm, O.D. 6-7µm, without spike (s)	002-5	20

Spermatid pipettes				
Identity colour: Blue	Product information	Ordering information	Article number	Box
	<b>Spermatid micropipettes</b> <ul style="list-style-type: none"> <li>- prepared from borosilicate glass tubing (O.D. 1mm, I.D. 0.78mm)</li> <li>- total length of the pipettes 5.5cm</li> <li>- bending angle 20°-45°</li> <li>- length of tip 0.5mm</li> <li>- bevelled 35°</li> <li>- with or without spike</li> <li>- with I.D. of the tip 7-8µm.</li> </ul>	I.D. 7-8µm, O.D. 9-10µm, bend 30 degrees, with spike (s)	002-7-30-s	20
		I.D. 7-8µm, O.D. 9-10µm, bend 30 degrees, without spike (s)	002-7-30	20
Hatching pipettes				
Identity colour: Pink	Product information	Ordering information	Article number	Box
	<b>Hatching micropipettes</b> <ul style="list-style-type: none"> <li>- prepared from borosilicate glass rod (O.D. 1.2mm),</li> <li>- total length of the pipettes 5.5cm</li> <li>- bending angle 20-45°</li> <li>- length of tip 1mm</li> <li>- with short thin taper S or with long thin taper L</li> <li>- sharp point</li> </ul>	with short taper, bend 35 degrees	003-35-S	20
		with long taper, bend 35 degrees	003-35-L	20
		with short taper, bend 30 degrees	003-30-S	20
		with long taper, bend 30 degrees	003-30-L	20
Biopsy pipettes				
Identity colour: Yellow	Product information	Ordering information	Article number	Box
	<b>Biopsy micropipettes</b> <ul style="list-style-type: none"> <li>- prepared from borosilicate glass tubing (O.D. 1mm, I.D. 0.78mm),</li> <li>- total length of the pipettes 5.5cm</li> <li>- bending angle 20°-45°</li> <li>- length of tip 1mm</li> <li>- with blunt opening A or bevelled 40° B</li> <li>- with I.D. 10,15, 20 and 35µm.</li> </ul>	I.D. 10µm, angle 35 degrees, blunt opening	004-10-35-A	20
		I.D. 10µm, angle 35 degrees, bevelled	004-10-35-B	20
		I.D. 15µm, angle 35 degrees, blunt opening	004-15-35-A	20
		I.D. 15µm, angle 35 degrees, bevelled	004-15-35-B	20
		I.D. 35µm, angle 35 degrees, blunt opening	004-35-35-A	20
		I.D. 35µm, angle 35 degrees, bevelled	004-35-35-B	20



## Advantages

- The pipettes are made of medical grade borosilicate glass.
- To meet international standards as well as the requirements of the FDA, the pipettes are sterilised by gamma-radiation.
- The pipettes are manufactured in accordance with the highest quality standards respecting the following parameters:
  - length of the pipette
  - bending angle of the tip
  - inner and outer diameter of the tip
  - length of the tip
- Upon customer request, we are prepared to produce any specification.
- By means of microbiological tests as well as by mouse-embryo-test (MEA) each lot is screened to assure the highest possible safety for the user.
- The pipettes are singly and sterile packed and can be taken out of the protective package easy, quick and safe.
- For better differentiation the pipettes are delivered in different coloured packages.

## Indices

- Production and packing of the pipettes happens under cleanroom conditions.
- Microtech produces in compliance with council directive 93/42/EEC.
- Gynemed distributes worldwide except Japan, Australia and Czech.
- For any Certificate of Analysis (i.e. MEA) look at [www.gynemed.de](http://www.gynemed.de)